

## Clemex ST 2.35 MP Colour USB 3.0 Camera **CL-13-875**

### **Sensor specification**

Image Sensor	Sony : IMX249 Colour
Shutter Type	Global Shutter
Optical Format	1/1.2"
Pixel Size	5.86 x 5.86 $\mu$ m
Resolution	1920 x 1200 pixels

### **Camera Specifications**

Frame Rate	40 fps at full resolution
Bit Depth	8 bit
Exposure Control	Software Control
Exposure Time	33 $\mu$ sec - 81 sec
Gain Range Analog	0 to 24 dB
Gain Range Digital	0 to 6 dB
White Balance	Manual and Auto control

### **Mechanical Specifications**

Data Interface	USB 3.0, Super Speed Micro B
Lens Mount	C-Mount
Dimensions	28 x 28 x 40 mm (excluding connectors)
Mass	50 g
Operating Temperature	0 to 50 °C
Storage Temperature	-30 to 70 °C
Operating Humidity	5 to 95 %, non-condensing
Shock / Vibration	50 G shock / 5 G (2 to 200Hz) vibration

### **Software, Power & Warranty**

Operating Systems	Windows 7, Windows 8, Windows 10 / 64-bit operating systems
Required Software	Clemex software version 8.0
Power Consumption	Less than 3 W
Power Requirement	5 V (Typ) (This conforms to the USB standard)
Warranty	One (1) year

